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百迈客生物科技  
BIOMARKER TECHNOLOGIES

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Species	Cells	Enzyme(s)	Reference
Human	脂肪细胞	2 型胶原酶: 0.01-0.5%	Effect of Collagenase Concentration on The Isolation of Small Adipocytes from Human Buccal Fat Pad., J Oral Sci , 2018
Human	脂肪细胞	Collagenase Type 1: 0.1%	Epigenome-wide Association Study of Body Mass Index, and the Adverse Outcomes of Adiposity, Nature541, 81, 2017
Human	基质细胞	胶原酶 2 型: 0.075%	Ultrasound-Assisted Liposuction Provides a Source for Functional Adipose-Derived Stromal Cells., Cytotherapy 19, 1491-1500, 2017
Human	骨髓基质细胞	胶原酶 1 型: 0.075%	Mesenchymal Stromal Cells Protect Human Cardiomyocytes from Amyloid Fibril Damage, Cytotherapy 19, 1426-1437, 2017
Human	间充质干细胞	胶原酶 2 型: 0.2%	Adipogenic Differentiation of Mesenchymal Stem Cells Alters Their Immunomodulatory Properties in a Tissue-Specific Manner., Stem Cells 35, 1636-1646, 2017
Human	脂肪源内皮细胞	胶原酶 1 型: 0.1%	Autologous Cell Sources in Therapeutic Vasculogenesis: In Vitro and In Vivo Comparison of Endothelial Colony-Forming Cells from Peripheral Blood and Endothelial Cells Isolated from Adipose Tissue., Cytotherapy 18, 242-52, 2016
Human	脂肪干细胞	胶原酶 1 型: 0.1%	High Glucose-Induced Reactive Oxygen Species Generation Promotes Stemness in Human Adipose-Derived Stem Cells, Cytotherapy 18, 371-83, 2016
Human	骨髓基质细胞	胶原酶 4 型: 0.2%	Effect of Mild Heat Stress on the Proliferative and Differentiative Ability of Human Mesenchymal Stromal Cells., Cytotherapy17,

			359-68, 2015
<b>Human</b>	脂肪干细胞	胶原酶 2 型: 0.1%	Expression Analysis of Human Adipose-Derived Stem Cells During In Vitro Differentiation to an Adipocyte Lineage., BMC Med Genomics 8, 41, 2015
<b>Human</b>	脂肪基质	胶原酶 1 型: 0.075%	Therapeutic Potential of Adipose-Derived SSEA-3-Positive Muse Cells for Treating Diabetic Skin Ulcers., Stem Cells Transl Med 4, 146, 2015
<b>Human</b>	脂肪组织细胞	胶原酶 1 型: 0.1%	Differential Effects of Processing Time and Duration of Collagenase Digestion on Human and Murine Fat Grafts., Plast Reconstr Surg 136, 189e-199e, 2015
<b>Human</b>	脂肪基质血管细胞	中性蛋白 酶: 2.4 u/ml、胶原 酶: 250 u /ml	Human White and Brite Adipogenesis is Supported by MSCA1 and is impaired by Immune Cells., Stem Cells 33, 1277-91, 2015
<b>Human</b>	脂肪源性间充质干细胞	胶原酶 2 型: 0.1%	Defined Serum-Free Media for In Vitro Expansion of Adipose-Derived Mesenchymal Stem Cells., Cytotherapy 16, 915, 2014
<b>Human</b>	脂肪提取干细胞	胶原酶 1 型: 0.15%	Choosing the Right Type of Serum for Different Applications of Human Adipose Tissue-Derived Stem Cells: Influence on Proliferation and Differentiation Abilities., Cytotherapy 16, 789, 2014
<b>Human</b>	间充质基质	胶原酶 1 型: 0.1%	Proliferative and Phenotypical Characteristics of Human Adipose Tissue-Derived Stem Cells: Comparison of Ficoll Gradient Centrifugation and Red Blood Cell Lysis Buffer Treatment Purification Methods., Cytotherapy 16, 1220-8, 2014
<b>Human</b>	脂肪间质干细胞	动物游离 胶原酶: 200u/ml	Xenofree Enzymatic Products for the Isolation of Human Adipose-Derived Stromal/Stem Cells., Tiss Eng 19, 473-8, 2013

<b>Human</b>	基质血管成分	胶原酶 1 型: 0.075%	Stromal Vascular Fraction Isolated from Lipo-Aspirates Using an Automated Processing System: Bench and Bed Analysis., J Tissue Eng Regen Med 7, 864, 2013
<b>Human</b>	脂肪干细胞	胶原酶 1 型: 0.1%	Platelet-Rich Plasma Greatly Potentiates Insulin-Induced Adipogenic Differentiation of Human Adipose-Derived Stem Cells Through a Serine/Threonine Kinase Akt-dependent Mechanism and Promotes Clinical Fat Graft Maintenance., Stem Cells Transl Med 1, 206-20, 2012
<b>Human</b>	血管周围干细胞	胶原酶 2 型: 0.1%	An Abundant Perivascular Source of Stem Cells for Bone Tissue Engineering., Stem Cells Transl Med 1, 673, 2012
<b>Human</b>	脂肪来源的间质血管	胶原酶 1 型: 0.1%	Concise Review: Adipose-Derived Stromal Vascular Fraction Cells and Platelet-Rich Plasma: Basic and Clinical Implications for TissueEngineering Therapies in Regenerative Surgery., Stem Cells Transl Med 1, 230-6, 2012
<b>Mouse</b>	脂肪基质	胶原酶 1 型: 0.1%	Adipose Stromal Vascular Fraction-Mediated Improvements at Late-Stage Disease in a Murine Model of Multiple Sclerosis., Stem Cells 35, 532-544, 2017
<b>Mouse</b>	脂肪间质细胞	胶原酶 2 型: 0.2%	Adipose Mesenchymal Stromal Cells Minimize and Repair Radiation-Induced Oral Mucositis.,Cytotherapy18, 1129-45, 2016
<b>Mouse</b>	脂肪基质	胶原酶 1 型: 0.1%	Improved Mobilization of Exogenous Mesenchymal Stem Cells to Bone for Fracture Healing and Sex Difference., Stem Cells 34, 2587-2600, 2016
<b>Mouse</b>	脂肪细胞	胶原酶: 0.1%	The Effects of A Single Developmentally-Entrained Pulse of Testosterone in Female Neonatal Mice On Reproductive and Metabolic Functions in Adult Life., Endocrinology 156, 3737, 2015

<b>Mouse</b>	间质血管 细胞	胶原酶 2 型: 0.2%	Natural Killer T Cells in Adipose Tissue are Activated in Lean Mice., Exp Anim 62, 319, 2013
<b>Mouse</b>	脂肪干细 胞	胶原酶 2 型: 0.1%	Biological and Clinical Availability of Adipose-Derived Stem Cells for Pelvic Dead Space Repair., Stem Cells Transl Med 1, 803, 2012

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- ❖ Dulbecco's Phosphate-Buffered Saline (DPBS)
- ❖ Hank's Balanced Salt Solution (HBSS)

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- ❖ Eagle's Minimum Essential Medium (EMEM) + 10% FBS
- ❖ Dulbecco's Modified Eagle Medium (DMEM) + 10% FBS
- ❖ Iscove's Modified Eagle Medium (IMEM) + 10% FBS
- ❖ Roswell Park Memorial Institute (RPMI) + 10% FBS
- ❖ Ham's F12 + 10% FBS
- ❖ 1:1 DMEM/F12 +10% FBS
- ❖ M199

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➤ Single-Cell RNA-Seq Reveals the Transcriptional Landscape and Heterogeneity of Aortic Macrophages in Murine Atherosclerosis

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➤ Developmental diversification of cortical inhibitory interneurons

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➤ Long-term correction of diabetes in mice by in vivo reprogramming of pancreatic ducts

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(BD Biosciences) 分选细胞。

➤ Graded Arrays of Spinal and Supraspinal V2a Interneuron Subtypes Underlie Forelimb and Hindlimb Motor Control

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➤ Conserved properties of dentate gyrus neurogenesis across postnatal development revealed by single-cell RNA sequencing

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➤ Single-cell profiling of human gliomas reveals macrophage ontogeny as a basis for regional differences in macrophage activation in the tumor microenvironment

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➤ Molecular Classification and Comparative Taxonomics of Foveal and Peripheral Cells in Primate Retina

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