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- [1] Yin P, Peter A, Franken H, et al. Preanalytical aspects and sample quality assessment in metabolomics studies of human blood[J]. Clinical chemistry, 2013: clinchem. 2012.199257.
- [2] Dunn W B, Broadhurst D, Begley P, et al. Procedures for large-scale metabolic profiling of serum and plasma using gas chromatography and liquid chromatography coupled to mass spectrometry[J]. Nature protocols, 2011, 6(7):

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- [1] Fernández-Peralbo M A, Castro L D. Preparation of urine samples prior to targeted or untargeted metabolomics mass-spectrometry analysis[J]. Trends in Analytical Chemistry, 2012, 41(41):75-85.
- [2] Bernini P, Bertini I, Luchinat C, et al. Standard operating procedures for pre-analytical handling of blood and urine for metabolomic studies and biobanks.[J]. Journal of Biomolecular Nmr, 2011, 49(3-4):231-243.
- [3] Yuille M, Illig T, Hveem K, et al. Laboratory management of samples in biobanks: European consensus expert group report.[J]. Biopreservation & Biobanking, 2010, 8(1):65.

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[1]Menezes E B, Velho A L C, Santos F et al. Uncovering sperm metabolome to discover biomarkers for bull fertility.[J]. BMC Genomics, 2019, 20: 714.

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[1]Bertoldo M J, Nadal-Desbarats L, Gérard N, et al. Differences in the metabolomic signatures of porcine follicular fluid collected from environments associated with good and poor oocyte quality [J]. Reproduction, 2013: REP-13-0142.

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[1] Junka A , Wojtowicz W , Z Bek A , et al. Metabolic profiles of exudates from chronic leg ulcerations[J]. Journal of Pharmaceutical and Biomedical Analysis, 2017, 137:13-22.

[2] Lam C W , Law C Y . Pleural Effusion Lipoproteins Measured by NMR Spectroscopy for Diagnosis of Exudative Pleural Effusions: A Novel Tool for Pore-Size Estimation[J]. Journal of Proteome Research, 2014, 13(9):4104-4112.

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[1]Wen X . Estradiol, progesterone, testosterone profiles in human follicular fluid and cultured granulosa cells from luteinized pre-ovulatory follicles[J]. *Reproductive Biology & Endocrinology* Rb & E, 2010, 8(1):117.2.8

[2]Bertoldo M J, Nadal-Desbarats L, Gérard N, et al. Differences in the metabolomic signatures of porcine follicular fluid collected from environments associated with good and poor oocyte quality[J]. *Reproduction*, 2013: REP-13-0142.

[3]Chao D L B J M , Boueilh T , Simard G , et al. Targeted metabolomics reveals reduced levels of polyunsaturated choline plasmalogens and a smaller dimethylarginine/arginine ratio in the follicular fluid of patients with a diminished ovarian reserve[J]. *Human Reproduction*, 2017, 32(11):2269-2278.

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- [1]Stoop, Marcel P., et al. "Quantitative proteomics and metabolomics analysis of normal human cerebrospinal fluid samples." *Molecular & Cellular Proteomics* (2010): mcp-M110.
- [2]Öhman A, Forsgren L. NMR metabonomics of cerebrospinal fluid distinguishes between Parkinson's disease and controls[J]. *Neuroscience letters*, 2015, 594: 36-39.
- [3] Tsutsui H, Muguruma Y, Noda T, et al. Development of targeted metabolomics for the determination of ornithine cycle compounds as possible biomarkers in cerebrospinal fluid regarding to Alzheimer's disease pathology using UHPLC-ESI-MS/MS[J]. *Med Mass Spectrom.*, 2018.

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on CF-associated viral and microbial communities.[J]. Journal of Cystic Fibrosis, 2013, 12(2):154-164.

[3] Garcia MR, Rodriguez JC, Navarro JF, et al. Molecular epidemiology of tuberculosis in elche, spain[J]. J Med Microbiol, 2002, 51(3): 273-274.

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[2]Sugimoto M ,Wong D T, Hirayama A, et al. Capillary electrophoresis mass spectrometry-based saliva metabolomics identified oral, breast and pancreatic cancer-specific profiles[J]. Metabolomics, 2010, 6(1):78-95.

[3]The human saliva metabolome[J]. Metabolomics, 2015, 11(6):1864-1883.

[4]Navazesh M . Methods for collecting saliva.[J]. Annals of the New York Academy of Sciences, 2010, 694(1):72-77.

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Untargeted Human Milk Metabolomics[J]. *Metabolites*, 2020, 10(2). [2]. Saben J L , Sims C R , Piccolo B D , et al. Maternal adiposity alters the human milk metabolome: associations between nonglucose monosaccharides and infant adiposity[J]. *The American Journal of Clinical Nutrition*. [3]. Sun H Z , Wang D M , Wang B , et al. Metabolomics of four biofluids from dairy cows: potential biomarkers for milk production and quality.[J]. *Journal of Proteome Research*, 2015, 14(2):1287-1298.

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[1]Nelson T M , Borgogna J C , Michalek R D , et al. Cigarette smoking is associated with an altered vaginal tract metabolomic profile[J]. *Scientific Reports*, 2018, 8(1):852.

[2]Kim TK, Thomas SM, Ho M, Sharma S, Reich CI, et al. (2009) Heterogeneity of vaginal microbial communities within individuals. *J Clin Microbiol* 47: 1181–1189

[1]NelsonTM,BorgognaJC,MichalekRD,etal.Cigarette smoking is associated with an altered vaginal tractmetabolomicprofile[J].*ScientificReports*,2018,8(1):852.

2.2.KimTK,ThomasSM,HoM,SharmaS,ReichCI,etal.(2009)Heterogeneityofvaginalmicrobialcommunitieswithinindividuals.*JClinMicrobiol*47:1181–1189

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[1] Dehui Xu, Yujing Xu et al. Alteration of metabolite profiling by cold atmospheric plasma treatment in human myeloma cells. *Cancer Cell Int.* 2018.

[2] Zuzana Racovaa, Eva Anzenbacherova, et al. Metabolite profiling of natural substances in human: in vitro study from fecal bacteria to colon carcinoma cells (Caco-2). *Journal of Nutritional Biochemistry.* 2020.

[3] Sarah Hayton, Arth L. Maker et al. Experimental design and reporting standards for metabolomics studies of mammalian cell lines. *Cellular and Molecular Life*

Sciences. 2017.

[4]Katja Dettmer,Nadine Nürnberger et al. Metabolite extraction from adherently growing mammalian cells for metabolomics studies: optimization of harvesting and extraction protocols.

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[1] Silas G.Villas-Boas, Extracellular metabolomics: A metabolic footprinting approach to assess Wber degradation in complex media.Analyticalbiotechnology, 2005.

[2] Hanna Meyer, Hendrikje Weidmann, Methodological approaches to help unravel the intracellular metabolome of Bacillus subtilis. Microbial Cell Factories2013, 12:69.

[3] Farhana R. Pinu, Analysis of Intracellular Metabolites from Microorganisms:uenching and Extraction Protocols, metabolites, 2017.

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[1]De Vos R C, Moco S, Lommen A, et al. Untargeted large-scale plant metabolomics using liquid chromatography coupled to mass spectrometry.[J]. Nature Protocols, 2007, 2(4):778-791.

[2]Kim H K, Choi Y H, Verpoorte R. NMR-based metabolomic analysis of plants[J]. Nature Protocols.

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- [1] Yuting Wang, Wenjie Ren, Yan Li, et al. Nontargeted metabolomic analysis to unravel the impact of di(2-ethylhexyl) phthalate stress on root exudates of alfalfa (*Medicago sativa*). *Science of the Total Environment* 2018.
- [2] Khorassani R, Hettwer U, Ratzinger A, Steingrobe B, Karlovsky P, Claassen N. 2011. Citramalic acid and salicylic acid in sugar beet root exudates solubilize soil phosphorus. *BMC Plant Biology* 11:121.
- [3] Luo Q, Wang S, Sun L N, et al. Metabolic profiling of root exudates from two ecotypes of *Sedum alfredii* treated with Pb based on GC-MS[J]. *Scientific Reports*, 2017, 7:39878.
- [4] Zhao L, Huang Y, Hu J, et al. ¹H NMR and GC-MS based metabolomics reveal defense and detoxification mechanism of cucumber plant under nano-Cu stress[J]. *Environmental Science & Technology*, 2016:acs.est.5b05011.

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[1] Bertini I, Hu X, Luchinat C. Global metabolomics characterization of bacteria: pre-analytical treatments and profiling [J]. *Metabolomics*, 2014, 10(2):241-249.

[2] Hanna Meyer, Hendrikje Weidmann, Methodological approaches to help unravel the intracellular metabolome of *Bacillus subtilis*. *Microbial Cell Factories* 2013, 12:69.

[3] Kathleen F Smart, Analytical platform for metabolome analysis of microbial cells using methyl chloroformate derivatization followed by gas chromatography– mass spectrometry, *nature protocols* | VOL.5 NO.10 | 2010.

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[1] The Preparation of Soil Samples for Chemical Analysis.G. M. MacNider; Journal of Industrial & Engineering Chemistry 1909 1 (7).

[2] Bartlett, R. J. and B. R. James. 1980. Studying dried, stored soil samples--some pitfalls. Soil Sci. Soc. Am. J.

[3] Stefan Jenkins, Tami L. Swenson. 2017. Construction of viable soil defined media using quantitative metabolomics analysis of soil metabolites. Frontiers in microbiology, 22.

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[1]David J. Beale .Microbial Metabolomics: Applications in Clinical, Environmental, and Industrial Microbiology [M]. Springer International Publishing, 2016: 160-161.

[2]Jun Y , Jun Z , Tao W , et al. Root exudates drive the soil-borne legacy of aboveground pathogen infection[J]. Microbiome, 2018 Sep 12;6(1).

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[1] Bi, H., Cai, D., Zhang, R., Zhu, Y., Zhang, D., Qiao, L., & Liu, Y. (2019). Mass Spectrometry-Based Metabolomics Approach to Reveal Differential Compounds in Pufferfish Soups: Flavor, Nutrition, and Safety. Food Chemistry, 125261.